
## Stanbio ALT/GPT (Liqui-UV®) Procedure No. 2930

For the Kinetic Quantitative Determination of ALT/  
GPT in Serum for Manual and/or Automated  
Procedures

### Summary and Principle

The enzyme alanine aminotransferase is widely reported in a variety of tissue sources. The major source of ALT is of hepatic origin and has led to the application of ALT determinations to the study of hepatic diseases. Elevated serum levels are found in hepatitis, cirrhosis, and obstructive jaundice. Levels of ALT are only slightly elevated in patients following a myocardial infarction.<sup>1</sup>

UV methods for ALT determination were first developed by Wroblewski and LaDue in 1956.<sup>2</sup> The method was based on the oxidation of NADH by lactate dehydrogenase (LDH). In 1980, the International Federation of Clinical Chemistry recommended a reference procedure for the measurements of ALT based on the Wroblewski and LaDue procedure. The Stanbio ALT reagent is based on a modified formulation of the IFCC<sup>3</sup> and Bergmeyer<sup>4</sup> methods.

Alanine aminotransferase (ALT) catalyzes the transfer of the amino group from alanine to 2-oxoglutarate, to form glutamate and pyruvate. The pyruvate formed is then reduced to lactate in the presence of LDH, with simultaneous oxidation of reduced nicotinamide adenine dinucleotide (NADH). The rate of decrease in absorbance at 340 nm is directly proportional to ALT activity.

### Reagent

#### ALT Buffer (R1), Ref. No. 2931

Composition:

|                     |            |
|---------------------|------------|
| L-Alanine           | 500 mmol/L |
| LDH (rabbit muscle) | 1200 U/L   |
| Tris Buffer, pH 7.5 | 100 mmol/L |

#### ALT Enzyme (R2), Ref. No. 2932

Composition:

|                      |             |
|----------------------|-------------|
| 2 - Oxoglutarate     | 15 mmol/L   |
| NADH (Disodium salt) | 0.18 mmol/L |

**Precautions:** The reagents are for "In Vitro Diagnostic Use". Normal precautions exercised in handling laboratory reagents should be followed. The reagents contain sodium azide which may be toxic if ingested. Sodium azide may also react with lead and copper plumbing to form highly explosive metal azides. Refer to Material Safety Data Sheet for any updated risk, hazard or safety information.

**Reagent Preparation:** Buffer and Enzyme liquid reagents are supplied ready-to-use. Prepare Working Reagent in the ratio of 5 parts Buffer (R1) to 1 part Enzyme (R2), (i.e., 25 mL Buffer and 5 mL Enzyme).

**Reagent Storage and Stability:** Reagents are stable until the expiration date on their respective labels, when properly stored at 2-8°C and protected from light. Reagents should appear clear and colorless. Discard if either appears cloudy or contains particulate matter. The Working Reagent is stable for 4 weeks at 2-8°C or 5 days at room temperature (15-30°C). The Working Reagent should be discarded if the initial absorbance, read against distilled water at 340 nm, is below 0.800.

### Material Required But Not Provided

Spectrophotometer capable of absorbance reading at 340 nm and 1 cm lightpath

Constant temperature block or bath, 37°C, or temperature controlled cuvette well

Accurate pipetting devices      Test tubes      Interval timer

### Specimen Collection and Storage

Non-hemolyzed serum is the specimen of choice. Whenever possible specimens should be separated and analyzed on the day of collection. Store serum in stoppered tubes. About 10% ALT is lost 3 days at 4°C and in 1 day at 25°C.<sup>5</sup>

**Interfering Substances:** Hemolysis must be avoided as the concentration of ALT in red cells is roughly 5 times that of serum.<sup>2</sup> Bilirubin levels up to 40 mg/dL and triglyceride levels up to 2000 mg/dL show no interference in this test.<sup>8</sup> Certain drugs and other substances are also known to affect ALT values.<sup>6</sup>

### Manual Procedure

1. Prepare ALT Working Reagent according to instructions.
2. Zero spectrophotometer at 340 nm with distilled water.
3. For each sample and control, add 1.0 mL Working Reagent to cuvette or test tube and warm to 37°C for 3 minutes.
4. Add 100 µL (0.10 mL) serum to its respective tube and mix gently.
5. Read and record absorbance at 1 minute. Continue incubating at 37°C and record absorbance again at 2 and 3 minutes. Rate should be constant.
6. Determine the average absorbance per minute (ΔA/min), multiply by factor -1746 for results in U/L.

NOTE: If cuvette is not temperature controlled, incubate samples at 37°C between readings.

**Quality Control:** Stanbio Ser-T-Fy I, Normal Control Serum, Cat. No. G427-86 and Stanbio Ser-T-Fy II, Abnormal Control Serum, Cat. No. G428-86 are recommended for each run. Other commercially available controls with ALT values assayed by this method are also suitable. ALT activity determined in these materials, by this procedure should fall within the ranges stated for the controls. Two levels of controls should be analyzed with each run.

**Calibration:** ALT activity is based on the "micromolar extinction coefficient" of NADH at 340 nm (see "Results" section). The instrument manufacturer's calibration guidelines should be followed to calibrate your analyzer. Assaying the ALT contents of a control serum with known ALT values can be used to assure instrument calibration has been performed correctly.

### Results

Values are derived based on the "absorptivity micromolar extinction coefficient" of NADH at 340 nm (0.0063). Units per liter (U/L) of ALT/GPT activity is that amount of enzyme which oxidizes one µmol/L of NADH per minute.

$$\text{U/L} = \frac{\Delta A / \text{Min}}{\text{Absorptivity}} \times \frac{\text{Total Volume}}{\text{Sample Volume}}$$

$$\text{U/L} = \frac{\Delta A / \text{Min}}{0.0063} \times \frac{1.10}{0.10}$$

$$\text{U/L} = \Delta A / \text{Min} \times 1746$$

### Limitations

If the ΔA/min. is greater than 0.342, dilute 1 part sample with 9 parts isotonic saline and re-assay. Multiply the result by 10. ALT values for neonatal patients have not been established with this procedure. Grossly icteric or turbid specimen may require the use of a sample blank.

### Expected Values<sup>8</sup>

Normal Range: 3 - 35 U/L (37°C)

This range should serve only as a guideline. It is recommended that each laboratory establish its own range of expected values, since differences exist between instruments, laboratories, and local populations.

### Performance Characteristics

**Comparison:** A group of 82 sera ranging in ALT activity from 3 - 682 U/L was assayed by the described ALT method and by a similar commercially available ALT reagent. Comparison of the results yielded a correlation coefficient of 0.999 and the regression equation was  $y = 0.988x + 0.43$ . (Comparison studies were performed according to NCCLS Tentative Guideline, EP9-T.)

**Precision:** Within-run precision was established by 20 assays on three different levels of commercial serum controls. Total Precision values were obtained by assaying the 3 commercial controls for 5 consecutive days.

|                      | Within-Run |         |         |
|----------------------|------------|---------|---------|
|                      | Serum 1    | Serum 2 | Serum 3 |
| Mean ALT (U/L)       | 22         | 45      | 102     |
| Std. Deviation (U/L) | 1.4        | 1.2     | 1.0     |
| C.V. (%)             | 6.2        | 2.6     | 1.0     |

|                      | Total Precision |         |         |
|----------------------|-----------------|---------|---------|
|                      | Serum 1         | Serum 2 | Serum 3 |
| Mean ALT (U/L)       | 23              | 43      | 99      |
| Std. Deviation (U/L) | 0.7             | 0.7     | 0.9     |
| C.V. (%)             | 3.1             | 1.6     | 0.9     |

Precision studies were performed according to NCCLS Tentative Guideline, EP5-T.

**Linearity:** Linear to 600 U/L at 37°C.<sup>8</sup> Performed according to NCCLS Guideline EP6-P.

**Sensitivity:** Based on an instrument resolution of A = 0.001, the method presented shows a sensitivity of 2.0 U/L.

### References

1. Henry, J.B.: Clinical Diagnosis and Management by Laboratory Methods, W.B. Saunders and Co., Philadelphia, PA, p-332-335 (1974)
2. Wroblewski, F. and LaDue, J.S. Proc. Soc. Exper. Biol. and Med. 91:569 (1956)
3. International Federation of Clinical Chemistry. Provisional Concentrations of Enzymes. Clin Chem 23: 887, 1977.
4. Bergmeyer HU, Scheibe P, Wahlefeld AW: Optimization of methods for aspartate aminotransferase and alanine aminotransferase. Clin Chem 24: 58, 1978.
5. Bergmeyer, H.U. Principles of Enzymatic Analysis. Verlag Chemic, 1978.
6. Young D.S. Effects of drugs on clinical laboratory tests. AACC Press, Washington D.C., 1990
7. Henry JB. Clinical Diagnosis and Management by Laboratory Methods, 17th ed. WB Saunders Co., 1984, p 1437.
8. Stanbio Laboratory Data

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### Referencias

1. Henry, J.B.: Clinical Diagnosis and Management by Laboratory Methods. W.B. Saunders and Co., Philadelphia, PA, p-332-335 (1974)
  2. Wroblewski, F. and LaDue, J.S. Proc. Soc. Exper. Biol. and Med. 91:569 (1956)
  3. International Federation of Clinical Chemistry. Provisional Concentrations of Enzymes. Clin Chem 23: 887, 1977.
  4. Bergmeyer HU, Scheibe P, Wahlefeld AW: Optimization of methods for aspartate aminotransferase and alanine aminotransferase. Clin Chem 24: 58, 1978.
  5. Bergmeyer, H.U. Principles of Enzymatic Analysis. Verlag Chemic, 1978.
  6. Young D.S. Effects of drugs on clinical laboratory tests. AACC Press, Washington D.C., 1990
  7. Henry JB. Clinical Diagnosis and Management by Laboratory Methods, 17th ed. WB Saunders Co., 1984, p 1437.
  8. Stanbio Laboratory Data
- Para asistencia técnica llame: 800-531-5535 • (830) 249-0772  
Fax (830) 249-0851 • e-mail: stanbio@stanbio.com  
http://www.stanbio.com  
Stanbio Laboratory • 1261 North Main Street • Boerne, Texas 78006  
DN: RBR.2930CE.01 • Last Revision: 01/04 • Procedure No. 2930CE