Stanbio Albumin LiquiColor® Procedure No. 0285

Quantitative Colorimetric Determination of
Albumin in Serum or Plasma

Summary and Principle

In 1964 bromocresol green (BCG) was reported as useful in the quantitative determination of serum albumin by Bartholomew and Delaney¹ and by Rodkey². This dye binding technique was again published a year later by Rodkey³ who described a BCG inverse colorimetric procedure specific for albumin at pH 7.05 and 615 nm. At the same time Watson⁴ reported BCG to be more sensitive for albumin than were either methyl orange or HABA (2-(4'-hydroxyazobenzene)-benzoic acid), two compounds also used in serum albumin quantitation by dye binding.

There followed reports by Dow and Pinto¹ and by Miyada et al⁶ of BCG method linearity up to 5 g/dL, minimum interference by lipids, hemoglobin and bilirubin, further evidence of high specificity and sensitivity, as well as excellent correlation for albumin values between BCG techniques and electrophoresis.

The method presented is essentially that of Dumas, Watson and Biggs^{7,8} modified by use of citrate instead of succinate buffer, a lower buffer concentration and reading of the final color at 550 instead of 630 nm⁹.

Reagents

Albumin Reagent, Cat. No. 0286

Contains bromocresol green, 18.8 mg/dL in citrate buffer, pH 4.2 ± 0.05 at 25°C. Also contains a surfactant.

Albumin Standard, Cat. No. 0287 (6.0 g/dL)

Aqueous solution of bovine fraction V albumin, equivalent to value stated on label. Sodium azide added as a preservative.

Precautions: For in Vitro Diagnostic Use.

Reagent Preparation: Reagent and Standard are ready to use.

Reagent Storage and Stability: Reagent is stable, if unopened and stored at 2 - 30°C, until expiration date on label. Standard is stable until expiration date on label when stored at 2-8°C. Bring Reagent and Standard to room temperature before use.

Materials Required But Not Provided

Spectrophotometer capable of absorbance readings at 550 nm (550-630 nm).
Accurate pipetting devices.

Cuvets Test Tubes
Interval timer

Specimen Collection and Preparation⁸

Serum, carefully removed from clot, is specimen of choice. Plasma may be used from blood collected with either EDTA or heparin as anticoagulant.

Sample Stability: Albumin is reportedly stable in serum for 1 week at room temperature, one month at 2-8°C and indefinitely when frozen.

Interfering Substances:¹⁰ BCG methods are reportedly free from interference by low levels of lipids, hemoglobin, bilirubin (20 mg/dL) and salicylate (60 mg/dL). Gross hemolysis must be avoided. Highly lipemic samples can be corrected by use of a serum blank.

Automated Analyzer

Parameters:

Wavelength	550 nm
Reaction Type.....	Endpoint
Reaction Direction	Increasing
Reaction Temperature	37°C
Sample/Reagent Ratio	1:100
Equilibration Time	3 Seconds
Read Time	4 Seconds
Lag Time	30 Seconds
Blank Absorbance Limit	0.450
High Absorbance	2.000A
Standard.....	6.0 g/dL
Low Normal	3.8 g/dL
High Normal	5.1 g/dL
Linearity	7.0 g/dL

Above parameters should be employed in programming automated analyzers for Albumin. Consult your instrument manual for programming instructions. Specific programming applications for most automated analyzers are available from Stanbio Customer Service Department.

Manual Procedure

1. Pipet into cuvetts the following volumes (mL) and mix well:

	Reagent Blank (RB)	Standard (S)	Sample (U)
Reagent	1.0	1.0	1.0
Standard	-	0.01	-
Sample	-	-	0.01

NOTE: All volumes may be increased x2 if the instrument requires volumes greater than 1.0mL.

2. Mix contents of each cuvet. Color develops immediately.

3. Read S and U vs RB at 550 nm within 15 minutes.

Quality Control: Use of commercial control serum, or pooled serum previously assayed and divided into frozen aliquots is recommended for use with each series of assays.

Results

Values are derived by the following equations:

$$1. \text{ Serum Albumin (g/dL)} = \frac{A_u}{A_s} \times 6.0$$

Where Au and As are the absorbance values of unknown and standard, respectively, and 6.0 the concentration of the standard (g/dL).

When a serum blank is required (icteric or hemolyzed specimen), label another tube SB. Add 1.0 mL "normal" saline, 0.01 mL serum, mix by inversion, transfer to cuvet and read absorbance (Asb) vs distilled water at 550 nm. Use this value to correct that of the unknown as follows:

$$2. \text{ Serum Albumin (g/dL)} = \frac{A_u - A_{sb}}{A_s} \times 6.0$$

NOTE: Samples having Albumin values greater than 7.0 dL are diluted 2-fold (1 + 1) with normal saline (sodium chloride, 8.5 g/L), the assay repeated and results multiplied by the dilution factor of 2.

Expected Values⁸

Normal Range 3.8-5.1 g/dL

Performance Characteristics⁹

Precision: Replicate assays (21) on each of 6 serum pools (range 3.2-4.3 g/dL) over a 17-day period showed standard deviations (SD) of 0.07-0.12 g/dL.

Correlation: Comparison of results derived on 20 patient sera (range 3.2-4.3 g/dL) by the method described and by the Coulter Automated Direct (Non-Blanked) procedure resulted in a correlation coefficient (r) of .9738.

Linearity: When performed as directed this method is linear from 0 to 7 g/dL.

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