## Stanbio Total Calcium LiquiColor® (Arsenazo III) Procedure No. 0155

For the Quantitative Colorimetric determination of Total Calcium in serum.

### Summary and Principle

Calcium measurements are used in the diagnosis and treatment of parathyroid disease, a variety of bone diseases, chronic renal disease and tetany (intermittent muscular contractions or spasma).

Accurate and precise measurement of calcium in biological fluids has traditionally been difficult<sup>1</sup>. Atomic absorption spectrophotometry (AAS) and the spectrophotometric measurement of calcium-dye complexes such as o-cresolphthalein complexone (CPC) are most often used for measuring total calcium. Variations of the latter method have been applied to automated analyzers.

The Stanbio method for the determination of calcium presented here uses Arsenazo III (3,6-bis[(2-Arsonophenyl) Azo]-4,5-dihydroxy-2, 7-naphthalenedisulfonic acid)<sup>2</sup>; CAS registry number: 1668-00-4. Arsenazo III is chemically stable and has a very high affinity for calcium in a neutral pH range. In this assay system, the Arsenazo III forms a 1:1 violet Arsenazo III: calcium complex with an absorbance maximum at 650 nm<sup>3</sup>. The concentration of calcium is proportional to the absorbance of the violet colored Arsenazo III: calcium complex. The color produced by this complex is stable for at least 30 minutes at room temperature (15 - 30°C).

Magnesium does not interfere in this assay system. The affinity of Arsenazo III for magnesium is essentially zero at the pH at which the assay is performed<sup>4</sup>.

### Reagents

#### Calcium Color Reagent, Cat. No. 0156

Buffered solution, 100 mmol/L, with Arsenazo III, 0.15 mmol/L, non-reactive stabilizers and surfactant.

#### Calcium Standard (10 mg/dL), Cat. No. 0157

Calcium carbonate in dilute hydrochloric acid.

#### Precautions: For In Vitro Diagnostic Use.

Do not pipet Color Reagent by mouth. Do not ingest. Organic arsenic compounds have been classified as potentially carcinogenic: therefore, safe laboratory practices should be carefully observed. Components of the reagent may be irritating to the skin and mucous membranes; avoid contact. If contact occurs, wash with copious amounts of water.

A major source of difficulty with the assay of calcium is contaminated glassware employed in the performance of the test. Many detergents and water supplies contain calcium and incompletely rinsed containers used in the test will lead to inaccurate results.

**Reagent Preparation:** The reagent and standard are supplied ready-to-use.

**Reagent Storage and Stability:** The reagent and standard are stable until their respective expiration dates when stored at room temperature (15 - 30°C). Reagent absorbance greater than 0.500 at 650 nm, when read against water, indicates deterioration and should therefore not be used.

### Materials Required But Not Provided

Spectrophotometer capable of absorbance readings at 650 nm  
Accurate pipetting devices  
Cuvets Vortex type mixer Timer

### Specimen Collection and Preparation

Clear, unhemolyzed serum is the preferred specimen. The specimen should be collected following the guidelines of NCCLS document H4-H3.<sup>5</sup> Heparinized plasma must not be used. Plasma prepared using EDTA, oxalate, citrate, which function by removal of calcium, obviously must not be used.

**Sample Stability:** Calcium is stable in serum for up to 24 hours at room temperature, approximately 1 week at 2-8°C, and for up to 1 year frozen at -20°C.<sup>6</sup> Specimens showing a precipitate should not be used, as calcium may have combined with fatty acids or proteins.

**Interfering Substances:** Contaminated glassware is the greatest source of error. All glassware should be washed in warm 20% nitric acid, rinsed well with distilled or deionized water, dried and stored in a dust-free area. Young et al.<sup>7</sup> have published a comprehensive list of drugs and substances which may interfere with in vitro diagnostic assays, including that for serum calcium. Interference from lipemia is minimized because of the small amount of sample used.

### Automated Analyzer

#### Parameters:

Wavelength.....	650 nm
Reaction Type.....	Endpoint
Reaction Direction.....	Increasing
Reaction Temperature.....	37°C
Sample/Reagent Ratio.....	1:100
Incubation Time.....	60 second
Blank Absorbance Limit.....	0.500 A
High Absorbance.....	2.000 A
Standard.....	10 mg/dL
Low Normal.....	8.5 mg/dL
High Normal.....	10.5 mg/dL
Linearity.....	15 mg/dL

Parameters should be employed in programming automated analyzers for Total Calcium. Consult your instrument manual for programming instructions. For a specific instrument application contact Stanbio's Customer Service Department.

### Manual Procedure

- Pipet into appropriately labelled cuvettes the following volumes (mL) and mix well: **Incubate for at least 60 seconds.**

	Reagent Blank (RB)	Standard (S)	Unknown (U)
Color Reagent	1.0	1.0	1.0
Specimen	—	—	0.01
Standard, (10 mg/dL)	—	0.01	—

- Read absorbance (A) of S and U against the RB at 650 nm. Final color is stable for at least 30 minutes.

NOTE: Volumes may be doubled if the instrument requires volumes greater than 1.0 mL.

**Quality Control:** Serum controls are recommended to monitor the performance of the manual and automated assay procedures, providing a continued screening of the instrument, reagents and techniques. Stanbio Ser-T-FY® I (Normal) and Ser-T-Fy® II (Abnormal) controls are commercially available for this use.

### Results

Results are derived by the following equation:

$$\text{Serum Calcium (mg/dL)} = \frac{\text{Au}}{\text{As}} \times 10$$

where **Au** and **As** are absorbance values of **Unknown** and **Standard**, respectively, 10 is the concentration of the **Standard**.

**Example:** Absorbance values of a serum sample was 0.581 and that of the standard was 0.645. Therefore :

$$\text{Serum Calcium (mg/dL)} = \frac{0.581}{0.645} \times 10 = 9.0 \text{ mg/dL}$$

### Expected Values<sup>7</sup>

Normal Range: Serum 8.5 - 10.5 mg/dL (2.1 - 2.6 mmol/L)

### Performance Characteristics<sup>4</sup>

**Precision:** Within-run precision was established by 9 assays on three different levels of serum specimens. Total Precision values were obtained by assaying the 3 same serum specimens for 5 consecutive days:

	Within-Run		
	Serum 1	Serum 2	Serum 3
Mean	6.8	12.9	21.8
Std. Deviation	0.005	0.117	0.13
C.V. (%)	0.07	0.9	0.7

	Total Precision		
	Serum 1	Serum 2	Serum 3
Mean	6.8	12.9	21.8
Std. Deviation	0.08	0.10	0.18
C.V. (%)	1.2	0.7	0.8

**Correlation:** Employing a Calcium o-Cresolphthalein as a reference method, 72 serum specimens ranging from 6.0 mg/dL to 18.0 mg/dL were assayed. The correlation coefficient was 0.998 and the regression equation was  $y = 0.99984x + 0.043$

**Sensitivity:** Using a 1:100 sample to reagent ratio and reading at 650 nm, a 1 mg/dL calcium sample will produce a net absorbance of approximately 0.065.

**Linearity:** This procedure is linear up to 15 mg/dL.

### References

- Teitz, N.W., Textbook of Clinical Chemistry. W.B. Saunders Co. (1986), pp. 1342.
- Bauer, P.J., Anal. Biochem. 110, 61-72, 1981
- Michaylova, V., and Ilkova, P., Anal Chim Acta 53, 194-198, 1971
- Stanbio Data
- NCCLS Document H4-H3, 3rd Ed. (1991)
- Tietz, N.W., editor, Fundamentals of Clinical Chemistry, 3rd. Ed., Saunders Co., 1987.
- Young, D.S., Pestaner, L.S. and Gibberman, V., Clin. Chem. 21. 1975.
- Kaplan, L.A., Pesce, A.J., Clinical Chemistry: Theory, Analysis, and Correlation, C.V. Mosby Co., St. Louis, 1984.

STANBIO LABORATORY, LP DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES OF THE MERCHANTABILITY AND FITNESS PERTAINING TO THIS PRODUCT WHICH ARE NOT EXPRESSLY DETAILED IN THIS PACKAGING INFORMATION OR A WRITTEN AGREEMENT BETWEEN THE BUYER AND SELLER OF THIS PRODUCT.

STANBIO LABORATORY, LP MAINTAINS THAT THIS PRODUCT CONFORMS TO THE INFORMATION CONTAINED IN THIS INSERT. PURCHASER MUST DETERMINE THE SUITABILITY OF THE PRODUCT FOR ITS PARTICULAR USE. USE ONLY IN ACCORDANCE WITH LABELING INSTRUCTIONS.

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8. Hoffman, W.S., The Biochemistry of Clinical Medicine, 4th Ed., 548-613, Year Book Medical Publishers, Inc., 1970.