


## Stanbio Cholesterol LiquiColor® Procedure No. 1010

Quantitative-Enzymatic-Colorimetric Determination  
of Total and HDL Cholesterol in Serum or Plasma

### Summary and Principle

The enzymatic approach to cholesterol methodology was introduced in 1973 by Flegg<sup>1</sup> and Richmond<sup>2</sup>, using cholesterol oxidase of bacterial origin following chemical saponification of the cholesterol esters. Roeschlau<sup>3</sup> modified this technique and Allain et al.<sup>4</sup> published the first fully enzymatic assay, combining cholesterol oxidase and cholesterol esterase. The method presented is based on the Allain procedure and utilizes these enzymes in combination with the peroxidase/phenol-4-antipyrine reagent of Trinder.<sup>5</sup>

Cholesterol esterase (CE) hydrolyzes esters to free cholesterol and fatty acids. The free cholesterol so produced plus the preformed cholesterol are then oxidized in the presence of cholesterol oxidase (CO<sub>x</sub>) to cholest-4-en-3-one and hydrogen peroxide. A quinoneimine chromogen, with absorption maxima at 500 nm, is produced when phenol is oxidatively coupled with 4-aminophenazone in the presence of peroxidase (POD) with hydrogen peroxide. The intensity of the final red color is proportional to total cholesterol concentration. Lipid Clearing Factor (LCF): a mixture of special additives developed by Stanbio is integrated into the cholesterol reagent to help minimize interference due to lipemia.

**A technique for high density lipoprotein (HDL) cholesterol utilizing Stanbio HDL Precipitating Reagent, Cat. No. 0599, must be procured separately.**

### Reagents

#### Enzymatic Cholesterol Reagent, Ref. No. 1011

Reagent contains the following active ingredients at stated concentrations

4-Aminophenazone	0.25	mmol/L
Phenol	25.0	mmol/L
Peroxidase	> 5.0	U/mL
Cholesterol Esterase	> 0.15	U/mL
Cholesterol Oxidase	> 0.1	U/mL

Buffers and Stabilizers

#### Cholesterol Standard, Ref. No. 1012 (200 mg/dL)

Buffered aqueous solution of cholesterol with stabilizers, surfactants and preservative.

**Precautions:** *For In Vitro Diagnostic Use Only.* Dispose of reagents in accordance with local requirements.

**Reagent Preparation:** Reagent and Standard are ready to use.

**Reagent Storage and Stability:** Reagent and Standard are stable stored at 2-8°C until expiration date on label. Once Reagent is opened, contamination must be avoided. Bring Reagent and Standard to room temperature before use.

### Materials Required But Not Provided

Spectrophotometer capable of absorbance readings at 500 nm.

Accurate pipetting devices

Cuvets, Test Tubes, Centrifuge, Interval Timer, Mixer (Vortex type)

Constant temperature bath, or block, 37° (Optimal)

### Specimen Collection and Preparation<sup>6</sup>

Blood should be collected following a 12-hour fast. Specimen may be serum, or plasma collected with EDTA as anticoagulant. Avoid hemolysis.

**Sample Stability:** Both total cholesterol and HDL cholesterol are reportedly stable 4 days at 2-8°C. Extended stabilities at -20°C are 3 months for "Total" and 7-14 days for "HDL". Whenever possible, specimens should be separated and analyzed on the day of collection.

**Interfering Substances:** Anticoagulants such as fluoride and oxalate will result in false low values. The test is not influenced by hemoglobin values up to 200 mg/dL or by bilirubin levels up to 10 mg/dL. However, interference from grossly icteric and heavily hemolyzed specimens is correctable by use of a serum/plasma blank (refer to "Results" section).

### Automated Analyzers

#### Parameters:

Wavelength	500 nm
Reaction Type	Endpoint
Reaction Direction	Increasing
Reaction Temperature	37°C
Sample/Reagent Ratio	1:100
Equilibration Time	3 Seconds
Read Time	4 Seconds
Lag Time	300 Seconds
Blank Absorbance Limit	0.30A
High Absorbance	2.000A
Standard	200 mg/dL
Low Normal	120 mg/dL
High Normal	310 mg/dL
Linearity	750mg/dL

Above parameters should be employed in programming automated analyzers for Cholesterol. Consult your instrument manual for programming instructions. Specific programming applications for most automated analyzers are available from Stanbio Customer Service Department

### Test Performance

1. Pipet into cuvetts the following volumes (mL) and mix well:

	Reagent Blank (RB)	Standard (S)	Sample (U)
Reagent	1.0	1.0	1.0
Standard	—	0.01	—
Sample	—	—	0.01

2. Incubate all cuvetts at 37°C for 5 minutes or at room temperature for 10 minutes.

3. Read S and U vs. RB at 500 nm within 60 minutes.

**Quality Control:** Use of commercial control serum, or pooled serum previously assayed and divided into frozen aliquots, is recommended for use with each series of assays.

### Results

Values are derived by the following equations:

$$\text{Serum Total Cholesterol (mg/dL)} = \frac{A_u}{A_s} \times 200$$

Where Au and As are the absorbance values of unknown and standard, respectively, and 200 the concentration of the standard (mg/dL).

When a serum blank is required (icteric or hemolyzed specimen), label another tube SB (Step 1, "Procedure" section). Add 1.0 mL "normal" saline 0.01 mL serum, mix by inversion, transfer to cuvet and read absorbance (Asb) vs. distilled water at 500 nm. Use this value to correct that of the unknown as follows:

$$\text{Serum Total Cholesterol (mg/dL)} = \frac{A_u - A_{sb}}{A_s} \times 200$$

NOTE: Samples having cholesterol values greater than 750 mg/dL are diluted 3-fold (1+2) with normal saline (sodium chloride, 8.5 g/L), the assay repeated and results multiplied by the dilution factor of 3.

### Expected Values<sup>7</sup>

Recent data are presented showing "normal ranges" according to age for total and HDL Cholesterol, and as to the risk of coronary heart disease (CHD) for HDL. Cholesterol expressed as percent of total cholesterol. It is recommended that each laboratory establish its own range of expected values, since differences exist between instruments, laboratories, and local populations.

Total Cholesterol	(mg/dL)
Cord	45-100
Newborn	53-135
Infant	70-175
Child	120-200
Adolescent	120-210
Adult	140-310
Desirable for Adults	140-200

Risk of CHD (HDL as % of total)		
Risk	Male	Female
Dangerous	< 7	< 12
High	7-15	12-18
Average	15-25	18-27
Below Average	25-37	
Probable Protection	> 37	> 40

HDL Cholesterol		
Age	Male	Female
0-14	30-65	30-65
15-19	30-65	30-70
20-29	30-70	30-75
30-39	30-70	30-80
> 40	30-70	30-85
Mean	45	55
Blacks: Approx. 10 mg/dL higher		

### Performance Characteristics

**Reproducibility:** A study was performed on a control serum (mean = 128 mg/dL) and on a patient's serum (mean = 367mg/dL), which consisted of a series of 5 assays on each of 5 days. Coefficients of variation (CV) were within run 2.5% and 1.0% and between-run (day-to-day) 3.5% and 2.9%, respectively.

**Correlation:** Determination of cholesterol by the procedure described (y) and by the Boehringer-Manheim "BMC Cholesterol Monotest" (x) on 66 sera (range 125-550 mg/dL) showed a correlation coefficient (r) of .991 and a regression equation of y = 1.04x-9.3.

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