## Stanbio LDH Liqui-UV® Procedure No. 2940

For the Quantitative Determination of Lactate Dehydrogenase in Serum and EDTA Plasma.

### Summary and Principle

Elevations in serum lactate dehydrogenase (LDH) occur from myocardial infarction, liver disease, pernicious and megaloblastic anemias, pulmonary emboli, malignancies, and muscular dystrophy<sup>1</sup>. A combined analysis of lactate dehydrogenase and creatine kinase, particularly the isoenzymes of each, provides a definite diagnosis of acute myocardial infarction.<sup>2</sup>

The procedure presented is essentially the Buhl and Jackson<sup>3</sup> modification of Wacker<sup>4</sup> which optimizes reaction conditions.

LDH specifically catalyzes the oxidation of lactate to pyruvate with the subsequent reduction of NAD to NADH. The rate at which NADH forms is proportional to LDH activity. The method described determines NADH absorbance increase per minute at 340 nm.

### Reagent

#### LDH Buffer (Reagent 1), Ref. No. 2941

L-Lithium Lactate	100mmol/L
2-Methyl-2-Amino-1-Propanol pH 8.8	600mmol/L

#### LDH Enzyme (Reagent 2), Ref.No. 2942

NAD	6mmol/L
-----	---------

**Precautions:** The reagents are for "In Vitro Diagnostic Use". Normal precautions exercised in handling laboratory reagents should be followed. The reagents contain sodium azide which may be toxic if ingested. Sodium azide may also react with lead and copper plumbing to form highly explosive metal azides. Refer to Material Safety Data Sheet for any updated risk, hazard or safety information.

**Reagent Preparation:** Buffer and enzyme liquid reagents are supplied ready-to-use. Prepare Working Reagent in the ratio of 5 parts Buffer (R1) to 1 part Enzyme(R2) (i.e., 25 mL Buffer and 5 mL Enzyme).

**Reagent Storage and Stability:** Reagents are stable until the expiration date on their respective labels, when properly stored at 2-8°C and protected from light. R1 and R2 should appear clear and colorless. Discard if either appears cloudy or contains particulate matter. The Working Reagent is stable for 2 weeks at 2-8°C or 1 day at room temperature (15-30°C).

### Materials Required But Not Provided

Spectrophotometer capable of absorbance reading at 340 nm and 1 cm lightpath • Constant temperature block or bath, 37°C, or temperature controlled cuvet well • Accurate pipetting devices  
Test tubes • Interval timer

### Specimen Collection and Preparation

Clear, unhemolyzed serum is the specimen of choice. Plasma may be used if collected with EDTA.

**Sample Stability:** Serum LDH activity is reportedly stable at 2-8°C for 3 days. Frozen samples show decreased activity of isoenzymes LD-4 and LD-5 and, therefore of total LDH values.

**Interfering Substances:** Heparin, citrate and oxalate inhibit the assay. Icteric and lipemic serum require a serum blank. Young et al.<sup>5</sup> have reviewed the effects of drug and other conditions of LDH levels.

### Automated Analyzers

#### Parameters:

Wavelength .....	340 nm
Reaction Type .....	Kinetic
Reaction Direction .....	Increasing
Reaction Temperature .....	37°C
Sample/Reagent Ratio .....	1:20
Equilibration Time .....	60 Seconds
Read Time .....	60 Seconds
Blank Absorbance Limit .....	0.450A
High Absorbance Change/Min .....	0.120 ΔA/Min.
Factor .....	3376
Low Normal .....	80 U/L
High Normal .....	285 U/L
Linearity .....	800 U/L
Cuvet Lightpath .....	1.0 cm

Above parameters should be employed in programming automated analyzers for LDH. Consult your instrument manual for programming instructions.

### Manual Procedure

1. Prepare LDH working reagent according to instructions.
2. Zero spectrophotometer at 340 nm with distilled water.
3. For each sample and control, add 1.0 mL Working Reagent to cuvet or test tube and warm to 37°C for 3 minutes.
4. Add 50 μL (0.050 mL) serum to its respective tube and mix gently.
5. Read and record absorbance at 1 minute. Continue incubating at 37°C and record absorbance again at 2 and 3 minutes. Rate should be constant.
6. Determine the average absorbance per minute (ΔA/min), multiply by factor 3376 for results in U/L.

**NOTE:** If cuvet is not temperature controlled, incubate samples at 37°C between readings.

**Quality Control:** Stanbio Ser-T-Fy® I, Normal Control Serum, Cat. No. G427-86 and Stanbio Ser-T-Fy® II, Abnormal Control Serum, Cat. No. G428-86 are recommended for verifying accuracy and precision. Other commercially available controls with LDH values assayed by this method are also suitable. LDH activity determined in these materials, by this procedure should fall within the ranges stated for the controls. Two levels of controls should be analyzed each day of testing.

**Calibration:** LDH activity is based on the "micromolar extinction coefficient" of NADH at 340 nm (see "Results" section). The instrument manufacturer's calibration guidelines should be followed to calibrate your analyzer.

### Results

Values are derived based on the "absorptivity micromolar extinction coefficient" of NADH at 340 nm (0.00622). A unit per liter (U/L) of LDH activity is that amount of enzyme which produces one μmol/L of NADH per minute.

$$U/L = \frac{\Delta A/Min}{Absorptivity} \times \frac{Total Volume}{Sample Volume}$$

$$U/L = \frac{\Delta A/Min}{0.00622} \times \frac{1.050}{0.050}$$

$$U/L = \Delta A/Min \times 3376$$

### Expected Values<sup>6</sup>

Males	80-285 U/L (37°C)
Females	103-227 U/L (37°C)

This range should serve only as a guideline. It is recommended that each laboratory establish its own range of expected values, since differences exist between instruments, laboratories, and local populations.

### Performance Characteristics<sup>7</sup>

**Comparison:** A group of 70 sera ranging in LDH activity from 64 - 504 U/L was assayed by the described LDH method and by a similar commercially available LDH reagent. Comparison of the results yielded a correlation coefficient of 0.999 and the regression equation was  $y = 1.108x + 1.582$ . (Comparison studies were performed according to NCCLS Tentative Guideline, EP9-T).

**Linearity:** Linear to 800 U/L. Samples exceeding this value should be diluted 2-fold (1 + 1) with deionized water, the assay repeated and results multiplied by 2.

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### Referencias

1. Wilkinson JH. Principles and Practice for Diagnosis Enzymology. Year Book Medical Publishers, 1976
2. Kachmar JR: Enzymes. In Fundamentals of Clinical Chemistry, NW Tietz, Editor, Saunders, Philadelphia, 1976, p 674.
3. Sacks HJ, Lanchantin GF: An elevation of serum transaminase in the jaundice state. Am J Clin Pathol 33:97, 1960
4. Karmen A: A note on the spectrophotometric assay of glutamic-oxalacetic transaminase in human blood. J Clin Invest 34:131, 1955
5. Bergmeyer HU, Scheibe P, Wahlefeld AW: Optimization of methods for aspartate aminotransferase and alanine aminotransferase. Clin Chem 24: 58, 1978.
6. Expert Panel of Enzymes of the International Federation of Clinical Chemistry: Part 3. Revised IFCC method for aspartate aminotransferase. Clin Chem 24:720, 1978.
7. Demetriou JA et al. In Clinical Chemistry - Principles and Technics 2nd ed. RJ Henry et al. Eds. Harper & Row, Hagerstown MD, 1974, p 873.
8. International Federation of Clinical Chemistry. Provisional Recommendations on IFCC Methods for the Measurement of Catalytic Concentrations of Enzymes. Clin Chem 23: 887, 1977.
9. Young D.S. Effects of drugs on clinical laboratory tests. AACC Press, Washington D.C., 1990
10. Henry JB. Clinical Diagnosis and Management by Laboratory Methods, 17th ed. WB Saunders Co., 1984, p 1437.
11. Stanbio Laboratory Data

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