Stanbio Total Protein LiquiColor® Procedure No. 0250

Quantitative Colorimetric Determination of
Total Protein in Serum

Summary and Principle

Protein molecules contain a large number of peptide bonds. When treated with copper ions (Cu²⁺) in alkaline solution, a colored complex is formed between the copper and the carbonyl and imine groups of these peptides. As a similar reaction occurs with biuret (the simplest of such compounds formed by heating urea), the term "biuret reaction" was adopted.

The method described is based on the reports of Weischelbaum¹ and Gornal et al.² The violet color developed is proportional to the number of peptide bonds in the protein and is nearly independent of the relative concentration of albumin and globulin.³

Reagents

Total Protein LiquiColor® Reagent, Cat. No. 0251

Solution of copper sulfate pentahydrate, 0.3 g/dL, in aqueous sodium hydroxide, 0.8 g/dL. Also contains potassium iodide and potassium sodium tartrate.

Total Protein Standard, Cat. No. 0256 (10 g/dL)

Aqueous solution of Bovine Albumin, Fraction V, with sodium azide as preservative.

Precautions: *For In Vitro Diagnostic Use.*

Standard contains sodium azide as a preservative. May react with copper or lead plumbing to form explosive metal azide build up.

Reagent Preparations: The reagent and standard are ready to use.

Reagent Storage and Stability: Total Protein kit is stable at 2-8°C until the expiration date on the label. Total Protein Reagent is stable at 2-30°C. Store standard at 2-8°C.

Materials Required But Not Provided

Spectrophotometer capable of absorbance readings at 550 nm
Accurate pipetting devices
Test Tubes
Cuvets
Vortex mixer

Specimen Collection and Preparation⁴

Carefully separate serum from clot. Avoid hemolysis. Do not use plasma, as fibrinogen gives increased results.

Sample Stability: Total protein in serum is reportedly stable for 1 week at room temperature and for 30 days at 2-8°C.

Interfering Substances: Gross hemolysis and lipemia must be avoided. This method is **not** to be used for the determination of protein in urine or spinal fluid, due to low protein levels and high concentrations of interfering compounds.

Automated Analyzer

Parameters:

Wavelength 550 nm
Reaction Type Endpoint
Reaction Direction Increasing
Reaction Temperature 37°C
Sample/Reagent Ratio 1:100
Equilibration Time 4 Seconds
Read Time 3 Seconds
Lag Time 300 Seconds
Blank Absorbance Limit 0.200A
High Absorbance 1.200A
Standard 10 g/dL
Low Normal 6.6 g/dL
High Normal 8.3 g/dL

Above parameters should be employed in programming automated analyzers for Total Protein. Consult your instrument manual for programming instructions.

Manual Procedure

1. Pipet into appropriately labeled cuvetts the following volumes (mL) and mix well:

	UNKNOWN (U)	STANDARD (S)	REAGENT BLANK (RB)
Reagent	1.0	1.0	1.0
Sample	0.01	–	–
Standard	–	0.01	–

Note: For instruments requiring volumes greater than 1.0 mL, reagent and sample volumes may be doubled.

2. Allow all tubes to stand for at least 5 minutes.
3. Read S and U vs RB at 550 nm within 1 hour.

Quality Control: Control sera of known total protein content should be included with each set of unknowns.

Results

Values are derived by the following equation:

$$\text{Serum Total Protein (g/dL)} = \frac{A_u}{A_s} \times 10$$

Where Au and As are the absorbance values of UNKNOWN and STANDARD, respectively, and 10 the concentration of STANDARD (g/dL).

Expected Values⁴

Normal Range: 6.6 - 8.3 g/dL

Performance Characteristics⁵

1. Precision: Between run \pm 5.0%
2. Accuracy: Average deviation \pm 0.2 g/dL.
3. Specificity: Reaction equal with both albumin and globulin.
4. Linearity: Reaction obeys Beer-Lambert Law from 0-10 g/dL.

References

1. Weichselbaum TE. Am J Clin Pathol 7:40, 1946
2. Gornal AG, Bardawil CJ, David MM. J Biol Chem 177:751, 1949
3. Weissman N, Schoenbach EB, Armistead EB. J Biol Chem 187:153, 1950
4. Cannon DC et al. IN Clinical Chemistry – Principles and Techniques, 2nd Ed. RJ Henry et al., Eds. Harper & Row, Hagerstown, MD, 1974, pp 411-421
5. Stanbio Laboratory data

STANBIO LABORATORY, LP DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES OF THE MERCHANTABILITY AND FITNESS PERTAINING TO THIS PRODUCT WHICH ARE NOT EXPRESSLY DETAILED IN THIS PACKAGING INFORMATION OR A WRITTEN AGREEMENT BETWEEN THE BUYER AND SELLER OF THIS PRODUCT.

STANBIO LABORATORY, LP MAINTAINS THAT THIS PRODUCT CONFORMS TO THE INFORMATION CONTAINED IN THIS INSERT. PURCHASER MUST DETERMINE THE SUITABILITY OF THE PRODUCT FOR ITS PARTICULAR USE. USE ONLY IN ACCORDANCE WITH LABELING INSTRUCTIONS.

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- Gornal AG, Bardawil CJ, David MM. J. Biol. Chem. 177:751, 1949.