
## Stanbio Triglyceride LiquiColor® Procedure No. 2200

For the Quantitative Enzymatic-Colorimetric  
Determination of Triglycerides in Serum or Plasma

### Summary and Principle

Measurement of triglyceride levels, when performed in conjunction with other lipid assays, proves useful in the diagnosis of primary and secondary hyperlipoproteinemia. Triglyceride concentrations are also of interest in following the course of diabetes mellitus, nephrosis, biliary obstruction and various metabolic abnormalities resulting from endocrine disturbances.<sup>1, 2, 3</sup>

- Glycerol and fatty acids are first formed by lipase action on the triglycerides
- Glycerol is then phosphorylated by adenosine-5'-triphosphate (ATP) to produce glycerol-3-phosphate (G-3-P) and adenosine-5' -diphosphate (ADP) in a reaction catalyzed by glycerol kinase (GK):

- The G-3-P is oxidized by glycerylphosphate oxidase (GPO) producing dihydroxyacetone phosphate (DAP) and hydrogen peroxide:

- Peroxide reacts with a 4-aminoantipyrine and 4-chlorophenol under the catalytic influence of peroxidase (POD) to form quinoneimine:

Lipid Clearing Factor (**LCF**): a mixture of special additives developed by Stanbio is integrated into the triglyceride reagent to help minimize interference due to lipemia.

### Reagents

#### Triglyceride LiquiColor® Reagent, Cat. No. 2201

ATP 2.0 mM, Magnesium Salt 15.0 mM, 4-Aminoantipyrine 0.5 mM, 4-Chlorophenol 4mM, Glycerylphosphate Oxidase 1500 U/L, Sodium Azide 0.01% Lipase 4000 U/L, Glycerol Kinase 400 U/L, Peroxidase 2000 U/L, Good's Buffer 50 mM, pH 6.7 ± 0.1.

#### Triglyceride Standard, 200 mg/dL, Cat. No. 2103

Contains glycerol with surfactant to yield 200 mg/dL triglycerides as triolein. Sodium azide 0.01% added as a preservative.

#### Precautions: For In Vitro Diagnostic Use.

Reagent and standard contain sodium azide as a preservative. May react with copper or lead plumbing to form explosive metal azides. Upon disposal, flush with large amounts of water to prevent azide build up.

**Reagent Preparation:** Reagent and standard are supplied ready to use.

**Reagent Storage and Stability:** Triglycerides reagent is stable until expiration date on label when stored at 2-8°C and protected from light. Once reagent is opened, contamination must be avoided. Bring reagent and standard to room temperature before use.

### Materials Required But Not Provided

Spectrophotometer capable of absorbance readings at 500 nm (492-530 nm)  
Accurate pipetting devices • Heat block or water bath, 37°C  
Interval timer • Cuvets

### Specimen Collection and Preparation<sup>2</sup>

**Sample Stability:** Triglycerides are reportedly stable for at least 10 days at 2-8°C. Do not store samples at 15-30°C as phospholipids may hydrolyze, releasing free glycerol and falsely elevating triglyceride values.

**Interfering Substances:** Blood collecting devices containing glycerol (glycerin) cannot be used, such as those having stoppers so lubricated. Gross hemolysis or high bilirubin values, will produce falsely elevated values. A number of drugs and substances affect the determination of triglycerides.

Interference from gross icteric and heavily hemolyzed specimens is correctible by use of serum/plasma blank (refer to "Results" section)

### Automated Analyzer

#### Parameters:

Wavelength .....	500 nm
Reaction Type .....	Endpoint
Reaction Direction .....	Increasing
Reaction Temperature .....	37°C
Sample/Reagent Ratio .....	1:100
Equilibration Time .....	3 Seconds
Read Time .....	4 Seconds
Lag Time .....	300 Seconds
Blank Absorbance Limit .....	0.500
High Absorbance .....	2.000A
Standard .....	200 mg/dL
Low Normal .....	30 mg/dL
High Normal .....	150 mg/dL
Linearity .....	1000 mg/dL

Above parameters should be employed in programming automated analyzers for Triglycerides. Consult your instrument manual for programming instructions. Specific programming applications for most automated analyzers are available from Stanbio Customer Service Department.

### Manual Procedure

- Pipet into cuvetts the following volumes (mL) and mix well:

	Reagent Blank (RB)	Standard (S)	Sample (U)
Reagent	1.0	1.0	1.0
Standard	-	0.01	-
Sample	-	-	0.01

**NOTE:** Volumes may be increased 2-fold if the instrument requires volumes greater than 1.0mL.

- Incubate all cuvetts at 37°C for 5 minutes, or incubate at room temperature for 10 minutes.
- Read S and U vs RB at 500 nm within 60 minutes.

**Quality Control:** Commercial control sera of known triglyceride content (determined by this method) should be included with each set of Unknowns. Controls should be run following the same procedure as for unknown.

### Results

Values are derived by the following equations:

$$1. \text{ Serum Triglyceride (mg/dL)} = \frac{\text{Au}}{\text{As}} \times 200$$

where Au and As are the absorbance values of unknown and standard, respectively and 200 the concentration of the standard (mg/dL).

When a serum blank is required (icteric or hemolyzed specimen), label another tube SB (Step 1. "Procedure" section). Add 1.0 mL "normal" saline, 0.01 mL serum, mix by inversion, transfer to cuvet and read absorbance (Asb) vs distilled water at 500 nm. Use this value to correct that of the unknown as follows:

$$2. \text{ Serum Triglycerides (mg/dL)} = \frac{\text{Au} - \text{Asb}}{\text{As}} \times 200.$$

NOTE: Samples having triglyceride values greater than 1000 mg/dL are diluted 5-fold (1 + 4) with normal saline (sodium chloride, 8.5 g/L), the assay repeated and results multiplied by the dilution factor of 5.

### Expected Values<sup>4</sup>

30 – 150 mg/dL

It is recommended that each laboratory establish its own normal range of expected values, since differences exist between instruments, laboratories, and local populations.

### Performance Characteristics<sup>5</sup>

**Reproducibility:** A study was performed on a control serum (mean = 57 mg/dL) and a patient pool (mean = 327 mg/dL) which entailed 10 determinations on each for 5 successive days. Coefficients of variation (CV) were within run 1.18% and 0.88% and between runs 1.96% and 0.88%, respectively.

**Correlation:** Determination of triglyceride levels by the procedure described (y) and by the GPO method of Boehringer Mannheim (x) on 57 sera (range 47 to 950 mg/dL) showed a correlation coefficient (r) of .999 and a regression equation of  $y = 1.029x - 7.46$ .

**Linearity:** When performed as directed the method is linear from 0 to 1000 mg/dL.

### References

- Fredrickson DS, Levy RI, Lees RS: New Engl J Med 276:34, 1967
- Dryer RL: IN Fundamentals of Clinical Chemistry. NW Teitz, Ed. WB Saunders, Philadelphia, 1970, p329.
- Wahlefeld, AW: IN Methods of Enzymatic Analysis, Vol 5, HU Bergmeyer, Ed. Academic Press, New York, 1974, pp 1831-1835.
- Scheletter G, Nussel E: Arbeitsmed Sozialmed Pracentimed 10:25, 1975
- Stanbio Laboratory data

STANBIO LABORATORY, LP DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES OF THE MERCHANTABILITY AND FITNESS PERTAINING TO THIS PRODUCT WHICH ARE NOT EXPRESSLY DETAILED IN THIS PACKAGING INFORMATION OR A WRITTEN AGREEMENT BETWEEN THE BUYER AND SELLER OF THIS PRODUCT.

STANBIO LABORATORY, LP MAINTAINS THAT THIS PRODUCT CONFORMS TO THE INFORMATION CONTAINED IN THIS INSERT. PURCHASER MUST DETERMINE THE SUITABILITY OF THE PRODUCT FOR ITS PARTICULAR USE. USE ONLY IN ACCORDANCE WITH LABELING INSTRUCTIONS.

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