Stanbio Total & Direct Bilirubin Procedure No. 0240, 0245

For the Quantitative Colorimetric Determination of Total and Direct Bilirubin in Serum or Plasma.

Summary and Principle

Most routine clinical procedures for the determination of serum bilirubin are based on the classical Diazo reaction of Ehrlich¹, which was first applied to the estimation of this bile pigment in serum in 1913 by Van den Berg and Snapper². An excellent review of many bilirubin methodologies is given by Winkelman and his co-workers³.

When a solution of sulfanilic acid in dilute hydrochloric acid is combined with one of sodium nitrite, nitrous acid is formed. This unstable acid, which must be freshly prepared for each group of assays, react to form diazotized sulfanilic acid.. The latter (also termed "p-benzenediazonium sulfonate") couples with bilirubin to produce azobilirubin, the color intensity of which is proportional to bilirubin concentration.

Since bilirubin standards suffer from a high degree of instability, the method presented utilizes a aqueous solution of N-1-naphthylethylenediamine dihydrochloride as a calibrator, according to Bilissis and Speer⁴. This pure stable aromatic amine couples with diazotized sulfanilic acid under the test conditions to give an azo dye with spectral qualities similar to azobilirubin. The calibrator supplied has been prepared to yield, for the technique described, a color equivalent to the labeled bilirubin level at 540 nm.

Reagents

Total Bilirubin Reagent, Cat. No. 0241

Sulfanilic acid, 32 mmol/L in dilute hydrochloric acid. Accelerator and stabilizer added.

Bilirubin Oxidant, Cat. No. 0262

Sodium nitrite, 20 mmol/L, aqueous. Stabilizer added.

Bilirubin Calibrator (10 mg/dL), Cat. No. 0263

N-1-naphthylethylenediamine dihydrochloride, 0.346 mmol/L.

Stabilizer added. Equivalent to 10 mg/dL Bilirubin in method described.

For Determination of Direct Bilirubin (Supplied Separately)

Direct Bilirubin Reagent, Cat. No. 0245

Sulfanilic acid, 32 mmol/L in dilute hydrochloric acid and stabilizer added.

Precautions: For In Vitro Diagnostic Use.

Reagents are caustic. Do not pipet by mouth.

Reagent Storage and Stability: Total Bilirubin reagent, Direct Bilirubin reagent, oxidant, and calibrator are stable when stored at 15-30°C until the expiration date on their respective labels. An eventual color change of Total Bilirubin reagent may occur but this color will not interfere with the assay. Discard oxidant if a dark yellow color develops.

Working Reagent: Add 1 drop of oxidant to every 1 mL of Total and/or Direct Bilirubin reagent. Invert gently 3-4 times to mix well. Working reagent is stable for 7 days when stored at 15-30°C.

Materials Required But Not Provided

Spectrophotometer capable of absorbance readings at 540 nm
Accurate pipetting devices, Cuvettes, Vortex mixer, Interval timer

Specimen Collection and Preparation^{3,5}

1. Serum, or plasma collected by use of any of the common anticoagulants, should be free of hemolysis and obtained in the fasting state to avoid chylous specimens.
2. Samples must be protected from both sunlight and white artificial light, as bilirubin is highly photolabile. Reportedly, about 50% of bilirubin can be lost when exposed to sunlight for 1 hour.

Sample Stability: Bilirubin is stable in serum or plasma 4-7 days at 2-8°C and for 3 months when frozen (-20°C).

Interfering Substances: Gross hemolysis will cause falsely high bilirubin values, due to the inhibition of the diazo reaction by oxyhemoglobin^{3,6}. An extensive list of drugs or other agents interfering with diazo bilirubin methodologies has been reported by Young et al⁷.

Turbid or lipemic serum may cause a falsely elevated bilirubin values. It is suggested to run a patient serum blank with each patient test.

Automated Analyzers

Parameters:	Total	Direct
Wavelength	540 nm	540 nm
Reaction Type	Endpoint	Endpoint
Reaction Direction	Increasing	Increasing
Reaction Temperature	37°C	37°C
Sample/Reagent Ratio	1:20	1:10
Equilibrium Time	3 seconds	3 seconds
Read Time	4 seconds	4 seconds
Lag Time	300 seconds	180 seconds
Blank Absorbance Limit	0.100A	0.100A
High Absorbance	2.000A	2.000A
Calibrators	10 mg/dL	10 mg/dL
Low Normal	0.0 mg/dL	0.0 mg/dL
High Normal	1.2 mg/dL	0.5 mg/dL
Linearity	20.0 mg/dL	10.0 mg/dL

Above parameters should be employed in programming automated analyzers for Bilirubin. Consult your instrument manual for programming instructions. Specific programming applications for most automated analyzers are available from Stanbio Customer Service Department.

Manual Procedure (Total Bilirubin)

1. Pipet into cuvettes labeled Reagent Blank (RB), K (Calibrator), and S (Specimen) the following volumes (mL), except Oxidant, which is 1 drop, mixing well after each addition:

	TOTAL BILIRUBIN			
	RB	K	S	S
"Total" Reagent	1.0	1.0	-	1.0
Oxidant (Drop)	1	1	-	1
Calibrator (K)	-	0.05	-	-
Sample	-	-	-	0.05
Water	0.05	-	-	-

2. Allow cuvettes to incubate at room temperature for at **least** 5 minutes.
3. Read K and S vs. RB at 540 nm within 30 minutes.

Manual Procedure (Direct Bilirubin)

1. Pipet into cuvettes labeled RB (Reagent Blank), K (Calibrator), Specimen Blank (SB), and S (Specimen) the following volumes (mL), except Oxidant, which is 1 drop, mixing well after each addition:

	DIRECT BILIRUBIN			
	RB	K	SB	S
"Direct" Reagent	1.0	1.0	1.0	1.0
Oxidant (Drop)	1	1	-	1
Water	0.1	-	-	-
Calibrator (K)	-	0.1	-	-
Sample	-	-	0.1	0.1

2. Allow cuvettes to incubate at room temperature for **exactly** 3 minutes.
3. **Immediately** read the K and S-SB vs. RB at 540 nm.

Notes:

1. Timing of 3 minutes for "Direct" samples after addition of serum is critical, as absorbance readings will increase slowly due to presence of the "Indirect" fraction. Timing of calibrator is not critical, as reaction is complete in 3 minutes.
2. For newborns, usually only total bilirubin is requested. However, presence of the conjugated (direct) form occurs with sufficient frequency in newborns that the direct bilirubin test can be utilized.

Quality Control: Control sera, assayed by this method should be included with each run of unknowns.

Results

Values are derived from the following calculations:

$$\text{Total or Direct Bilirubin (mg/dL)} = \frac{\text{As-Asb}}{\text{Ak}} \times 10$$

Where As, Asb and Ak are the absorbance values of the Specimen, Specimen Blank and Calibrator respectively, and 10, the concentration of the calibrator.

Expected Values^{3,5}

Upper Limits of Normal (mg/dL):

Total Bilirubin	Adults:	1.2	Newborns:	12.0
Direct Bilirubin	Adults and Infants (over 1 month)	0.5		

Performance Characteristics⁸

Precision: Serum samples were assayed 20 times (within run precision) and for 10 working days (run-to-run precision). Coefficients of variation were "within run" 3.3% (mean = 0.49 mg/dL), 2.3% (mean = 10.2 mg/dL) and "run-to-run" 3.5% (mean = 0.51 mg/dL), 2.4% (mean = 10.2 mg/dL) for total bilirubin. "Within run" coefficient of variation was 1.82% (mean = 3.70 mg/dL) and 2.76% (mean = 3.62 mg/dL) for direct bilirubin (abnormal).

Correlation: Determination of bilirubin by the procedure described (y) and by a modified Jendrassik-Groff⁹ technique (x) on 25 sera showed a correlation coefficient (r) of 0.9995 for total (normal and abnormal) bilirubin and 0.9970 for direct bilirubin (abnormal). The regression equation for total bilirubin was y = 1.06x - 0.100.

Sensitivity: Both Total and Direct procedures showed a sensitivity of 0.018 mg/dL per 0.001 absorbance units.

Linearity: Total and Direct procedures have been established to 20 mg/dL and 10 mg/dL, respectively.

References

1. Ehrlich P. Charite Ann 8:140, 1883.
2. Van den Bergh AAH, Snapper J. Deut Arch Klin Med 110:540, 1913.
3. Winkelman JW et al. In Clinical Chemistry – Principles and Technics. RJ Henry et al., Eds. Harper & Row, Hagerstown, MD. 1974, pp1042-1061.
4. Bilissis PK, Speer RJ. Clin Chem 9:552, 1963.
5. Routh JI. In Fundamentals of Clinical Chemistry. NW Tietz, Ed Saunders, Philadelphia, 1976, pp. 1035-1043.
6. Shull DC. Clin Chem 26:22, 1980.
7. Young DS et al. Clin Chem 21:22, 1980.
8. Stanbio Laboratory data.
9. Jendrassik L., Groff, P. Biochem Z 297-81, 1938.

STANBIO LABORATORY, LP DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES OF THE MERCHANTABILITY AND FITNESS PERTAINING TO THIS PRODUCT WHICH ARE NOT EXPRESSLY DETAILED IN THIS PACKAGING INFORMATION OR A WRITTEN AGREEMENT BETWEEN THE BUYER AND SELLER OF THIS PRODUCT.

STANBIO LABORATORY, LP MAINTAINS THAT THIS PRODUCT CONFORMS TO THE INFORMATION CONTAINED IN THIS INSERT. PURCHASER MUST DETERMINE THE SUITABILITY OF THE PRODUCT FOR ITS PARTICULAR USE. USE ONLY IN ACCORDANCE WITH LABELING INSTRUCTIONS.

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