Stanbio Direct Creatinine LiquiColor® Procedure No. 0420

Quantitative Kinetic Determination of Creatinine
in Serum or Urine

Summary and Principle

In 1886 Jaffe¹ described a method for the determination of creatinine involving a protein free filtrate and reaction with picric acid in alkaline solution. Although several methods have been described since, the classic Jaffe reaction is still most widely used. This reaction is subject to interferences by a number of substances, including protein and glucose.^{2,3,4} Modifications have been developed to combat the drawbacks.⁵ Kinetic procedures have become popular because they are rapid, simple and avoid interference. The present method is based on a modification of the automated reaction rate of Fabinay and Eringshausen.⁶

Creatinine + Picric Acid Creatine-picric acid complex

Creatinine reacts with the picric acid in alkaline conditions to form a color complex which absorbs at 510 nm. The rate of color formation is proportional to the creatinine concentration in the sample.

Reagents

Creatinine Acid Reagent, Cat. No. 0421

Aqueous solution of picric acid, 3.6 mmol/L.

Creatinine Base Reagent, Cat. No. 0422

Aqueous solution of sodium hydroxide, 240 mmol/L.

Creatinine Standard - 5 mg/dL, Cat. No. 0423

Aqueous solution of creatinine zinc chloride.

Precautions: *For In Vitro Diagnostic Use Only.*

Do not ingest. Toxicity has not been established. Care should be exercised in handling the Creatinine Base Reagent since it is caustic. Picric Acid is a strong oxidizing agent. Avoid contact with

skin. Wipe away any spills, since dry picric acid is explosive.

Reagent Preparation: Mix 1 part Creatinine Acid reagent with 1 part Creatinine Base Reagent to prepare "Working Reagent". Standard is ready to use.

Reagent Storage and Stability: Creatinine Acid Reagent, Creatinine Base Reagent and Creatinine Standard are stable until the date of their expiration when properly stored at 15-30°C.

Working reagent is stable for 30 days when properly stored at 15-30°C. Do not use Working reagent if turbid.

Materials Required But Not Provided

Spectrophotometer capable of absorbance readings at 510 nm,
Constant temperature block or bath (37°C)

Temperature controlled cuvette well (37°C)

Interval timer

Accurate pipetting devices

Test Tubes

Vortex mixer

Specimen Collection and Preparation

Serum: Remove specimen from clot promptly to prevent hemolysis. Urine: Dilute a portion of well mixed and measured 24-hour urine collection 1:20 (1+19) with distilled water. Mix well.

Sample Stability: Creatinine values have a reported stability of one day at 2-8°C, and several months when frozen (-20°C) and protected from evaporation and contamination. Urine specimens must be preserved with 15 grams of boric acid. Urine is stable for 4 to 7 days at room temperature (15-25°C).

Interfering Substances: A number of substances affect the accuracy of creatinine determination. See Young, et al⁷ for a comprehensive list.

Automated Analyzers

Parameters:

Wavelength	510 nm
Reaction Type	Kinetic
Reaction Direction	Increasing
Reaction Temperature	37°C
Sample/Reagent Ratio	1:20
Lag Time	20 Seconds
Read Time	60 Seconds
Blank Absorbance Limit	<0.200
High Absorbance Change/Min	0.625ΔA/Min.
Standard Value	5 mg/dL
Low Normal	0.7 mg/dL
High Normal	1.5 mg/dL
Linearity	20 mg/dL

Above parameters should be employed in programming automated analyzers for Creatinine. Specific programming applications for most automated analyzers are available from Stanbio Customer Service Department.

Manual Procedure

1. Set spectrophotometer cuvette temperature to 37°C.
2. Zero spectrophotometer with D.I. water at 510 nm.
3. Pipet 1.0 mL of working reagent into test tubes and prewarm at 37°C for 3 minutes.
4. Transfer 0.050 mL (50 μL) of standard to test tube, mix and immediately place in cuvette well.
5. After exactly twenty (20) seconds read and record absorbance (A₁).
6. At exactly sixty (60) seconds after reading (A₁), read and record absorbance (A₂).
7. Calculate change in absorbance (ΔA) by subtracting (A₂-A₁).
8. Run patients and controls by following steps 4 through 7.

Quality Control: Two (2) levels of control material with known Creatinine levels determined by this method, should be analyzed each day of testing.

Results

Values are derived by comparing the absorbance change (ΔA) of the unknown (u) with that of a standard (s) identically treated.

$$\text{Serum Creatinine (mg/dL)} = \frac{\Delta A_u}{\Delta A_s} \times 5$$

Where ΔA_u and ΔA_s are the absorbance changes (increase) of unknown and standard, respectively, and 5 the concentration of standard (mg/dL).

Example: ΔA_u = 0.01, ΔA_s = 0.05

$$\text{Serum Creatinine (mg/dL)} = \frac{0.01}{0.05} \times 5 = 1.0$$

$$\text{Urine Creatinine (mg/dL)} = \frac{\Delta A_u}{\Delta A_s} \times 5 \times 20 \text{ (Dilution factor)}$$

$$\text{Urine Creatinine (mg/L)} = \text{Urine Creatinine (mg/dL)} \times 10$$

$$\text{Urine Creatinine (mg/L/24hrs.)} = \text{Urine Creatinine (mg/L)} \times \frac{24 \text{ hrs. (mL)}}{1000}$$

Expected Values⁸

Normal Range: Male (serum): 0.9 - 1.5 mg/dL
 Male (urine): 1000 - 2000 mg/24hrs.
 Female (serum): 0.7 - 1.4 mg/dL
 Female (urine): 600 - 1500 mg/24hrs.

This range should serve only as a guideline. It is recommended that each laboratory establish its own range of expected values, since differences exist between instruments, laboratories and local populations.

Performance Characteristics⁹

Correlation: Serum and control samples (n = 38) having creatinine values between 0.6 – 25 mg/dL were assayed by this method and by a similar commercially available kinetic procedure. Statistical analysis revealed a correlation coefficient (r) of 0.998, with a regression equation of y = 0.96x + 0.046.

Precision:	Mean (mg/dL)	SD	CV%
	Within Run		
	0.9	0.08	4.0
	8.6	0.09	2.1
	Between Run		
	1.1	0.08	4.8
	8.4	0.16	2.6

Sensitivity: Based on an instrument resolution of A = 0.001, the method presented shows a sensitivity of 0.16 mg/dL.

Linearity: Linear to 20 mg/dL. Samples exceeding this value should be diluted 2-fold (1+1) with deionized water, the assay repeated and results multiplied by 2.

References

1. Jaffe M. *Physiol Chem* 10:391, 1886
 2. DiGiorgio J. In *Clin. Chem. - Principles and Technics*, RJ Henry et al, Eds, Harper & Row, Hagerstown, MD, 1974, pp549 - 550.
 3. Cook JGH. *Ann Clin Biochem* 12:219, 1975
 4. Tausky HH. In *Standard Methods of Clinical Chemistry*, Vol 3, D Seligson, Ed. Academic Press, New York, 1966, p99
 5. Heinegard D, Tiderstrom G. *Clin Chim Acta* 43:305, 1973
 6. Fabinay DL, Eringshausen G. *Clin Chim* 17:696, 1971
 7. Young DS et al. *Clin Chem* 21:286 D, 1975 (Special Issue)
 8. Larsen K. *Clin Chim Acta* 41:209, 1972
 9. Stanbio Laboratory Data
- STANBIO LABORATORY, LP DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES OF THE MERCHANTABILITY AND FITNESS PERTAINING TO THIS PRODUCT WHICH ARE NOT EXPRESSLY DETAILED IN THIS PACKAGING INFORMATION OR A WRITTEN AGREEMENT BETWEEN THE BUYER AND SELLER OF THIS PRODUCT.
- STANBIO LABORATORY, LP MAINTAINS THAT THIS PRODUCT CONFORMS TO THE INFORMATION CONTAINED IN THIS INSERT. PURCHASER MUST DETERMINE THE SUITABILITY OF THE PRODUCT FOR ITS PARTICULAR USE. USE ONLY IN ACCORDANCE WITH LABELING INSTRUCTIONS.

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